



PRODUCT SHEET

Bariatric Stand Tall Walker

Strong and functional walker developed for the bariatric user.

Stand Tall Walker is intended for users weighing up to 325 kg.



- Lockable hand brakes
- Possible mounting of an IV-holder
- Brakes on rear wheels
- Padded forearm support
- Weight capacity up to 325 kg



Bariatric Stand Tall Walker

General information

This product sheet contains important information about the use of Stand Tall Walker and the target group for the product. It also includes a range of specifications on Stand Tall Walker.

The most recent version of this product sheet can be found on our website. If you need a product sheet written in a larger font size, the product

sheet is available in PDF format at cobirehab.com, which can be scaled up to the required font size.

You can always contact Cobi Rehab on + 45 7025 2522 or cobi@cobi.dk if you want more information about Stand Tall Walker.

Cobi Rehab reserves the right to change the product specifications without providing prior notice.

Use

Stand Tall Walker is used in hospitals and care centers where it provides support and safety while walking for users up to 325 kg. Stand Tall Walker is only suited for indoor use.

Stand Tall Walker should be used on a level, hard surface with a maximum slope of 3 degrees.

It is possible to mount an IV holder (or a urine-bag) on the Stand Tall Walker. The IV pole attachment has a length of 1010 mm and is attached to the outside of the walking frame.

The attachment of the IV pole increases the width of the Stand Tall Walker by max. 160 mm. The IV holder may only be mounted onto the Stand Tall Walker by authorized staff.

Construction

Stand Tall Walker consists of a wide, strong frame of 32 mm tubes and four massive wheels, of which the rear two wheels have brakes.

Stand Tall Walker has lockable hand brakes, padded forearm supports, and a side plate for mounting e.g., a urine bag.

The lockable hand brakes on the hand bar limit high speed in e.g., turns. The parking brakes can be locked when the user must stand up or sit down.

Stand Tall Walker is equipped with four horizontal, continuous bolts (two on each side) for height adjustment of forearm supports and hand rests.

The bolts are locked with nuts. The height of forearm supports, and hand rests can be adjusted from 1000-1300 mm without the use of tools.

The cup-shaped forearm padding can be easily adjusted in width for the best possible support and relief.

Inner handles

By using the inner handles, the user can slowly stand up from the edge of the bed or chair while maintaining the stability and control of the aid. Together, it reduces the risk of falls.

Cleaning

Stand Tall Walker can be cleaned with water and ordinary detergent.

Disinfection

Stand Tall Walker can be disinfected with ethanol 70-85% or a chlorine product 1000-1200 ppm.

Long term storage

Stand Tall Walker should be stored at temperatures between -10° and 50°C.

Specifications

Maximum user weight	325 kg
Total length	730 mm
Total width	800 mm
Total width (with IV-holder)	920 mm
Width between underarm support, from center to center, adjustable	590-690 mm
Turning radius	2 meters
Height, adjustable in 25 mm intervals	1000-1300 mm
Weight	20 kg
Massive wheel	Yes, 100 mm
Frame	Powder-lacquered chromium molybdenum steel
Underarm support	PU
ISO classification	ISO 12 06 12 - 06
Expected product service life	7 years
Risk class	I
UDI-DI	5740000102520
Grundlæggende UDI	5740000100009NC
The surface temperature of the Stand Tall Walker can increase or decrease when exposed to external sources of heat or cold.	

Ordering

Model	Item number
Bariatrisk Stand Tall Walker	0180-076-000

