



BRAKE THE ISOLATION

MOBIBED - a mobile bed

MobiBed is a flexible product which quickly can be adjusted to the user's need for social participation, rest, sitting or lying down.

The mobile rest trolley can increase the quality of life for many users. Safely and therapeutically, it allows the user to participate in social activities with an easy possibility for rest and safe positioning whenever needed. For staff it is easy to drive and maneuver. Angles and height can be adjusted with a remote control.



- ✓ Possible to adjust height, back and foot support with the remote control.
- ✓ The height ensures better working conditions for the staff.
- ✓ Used in institutions in Scandinavia since 2017.
- ✓ Cushions for individual positioning.
- ✓ Additional positioning cushions available.
- ✓ Three length settings in the back part and three in the seat part.
- ✓ Recommended up to 185 cm user height.
- ✓ Bedrails and belts available.
- ✓ Puncture free tires.
- ✓ Quiet, durable control system.
- ✓ Anti-tip.

MobiBed – individual positioning and mobility

Being part of the community can increase the quality of life

MobiBed is used to position people with very reduced mobility sitting or lying down (also lying on their side) creating possibilities to participate in activities and socializing.

MobiBed is equipped with four wheels with puncture-free tyres. The front wheels can turn 360°, which makes it easy to maneuver around. The rear wheels have brakes.

Bedrails are available in the same design as the frame. The bedrails increase the safety for the user and the control of the user for the care staff.

MobiBed's raising and lowering functions as well as adjusting the back angle, are activated using a remote control that can be attached to the control handle.

The maximum height of MobiBed measured from the floor to the hanger (in horizontal position) is 85 cm, which allows care staff to handle the user ergonomically correctly for the benefit of the working environment.



The MobiBed can be adjusted in the backrest and footrest, from a horizontal position to a 60° seated position.

There are several position options for the back part and the seat part.

Construction - MobiBed

The MobiBed is made of modules and cushion components that fit together without annoying seams.

MobiBed is intended for indoor use but can also be used outdoors on flat terrain such as terraces.

The recycling value is high, as the amount of cushion accessories and the flexibility of the modular structure make it easy to adapt the MobiBed to another user.

The maximum user weight is 90 kg and must be observed to avoid damage to the product.

MobiBed is not suitable for people who can move around on their own.



For more information on MobiBed and our other products

Please contact Cobi Rehab cobi@cobi.dk



COBI REHAB

| Fuglebækvej 1D, DK-2770 Kastrup, Denmark, +45 7025 2522, cobirehab.com