



# Bariatric wheelchair Transporter

Rigid steel frame wheelchair constructed to suit bariatric users

The bariatric wheelchair Transporter is a part of our XXL-Rehab series and is designed for easy indoor transport inside hospitals. The wheelchair can easily be disinfected, and it is free of maintenance.

The XXL-Rehab Transporter is designed for easy access and safety for the bariatric user and the caregiver. The 24" rear wheels are solid, provides stability and allows easy steering and movement of the chair using minimal energy.

The wheelchair has a basket attached under the seat for storing personal belongings during transit.

## The inside width

The inside width of the XXL-Rehab Transporter is 71 cm. Using a 5 cm high foam seat and 7 cm free space to the armrests allows space for adipose tissue.

## Solid footrests

The footrests can carry up to 100 kg each and are height and side adjustable. The side adjustable feature makes access to intimate hygiene possible and it ensures that bariatric users with adipose tissue on the inside of the thighs can sit in the wheelchair with legs apart (up to 110 cm) and still have foot support. The footrests are covered with foam.

## Functional armrests

The ability for bariatric users to reach far back to reach the armrest is usually not present. Therefore, the wheelchair is designed with extra-long armrests with a length of 64 cm that sticks out 7 cm from the seat base. This helps the bariatric user when getting into the chair. The swing away armrests can be removed.

## Backrest

The space between backrest and seat has a clearance of 32 cm. The clearance provides good support to the lumbar spine for users with a protruding behind. It also enables the caregivers to reach for a sling under the user when positioning the user.

- For indoor transport
- Easy to disinfect
- No maintenance
- Large 24" rear wheels with parking brakes
- Solid armrests

MAX 325 kg

Basic UDI DI  
5740000100009NC



## Bariatric wheelchair Transporter Specifications

Bariatric wheelchair Transporter	
Maximum weight capacity	325 kg
Seat width with standard armrests	71 cm
Total outer width with standard armrests	78 cm
Seat depth	56 cm
Seat height	55 cm
Seat to footrest	40-42-44-46-48 cm
Width of footrest (total width changes)	58-110 cm
Distance from floor to pushbar	100 cm
Width between armrests	54 cm
Height of backrest	48 cm
Height of footrests	40-48 cm
Distance from floor to seat underside	44.5 cm
Total length	129 cm
Total weight	30 kg
Front wheels	125 mm
Rear wheels	24" rear wheels with parking brakes
Maximum slope	3°
Seat	Pressure moulded PUR foam
Fire retardant seat	Yes
Frame	Powder-lacquered chromium molybdenum steel
Footrest cover	PU foam
ISO - classification	ISO 09 12 03
Basic UDI	5740000100009NC
Risc class	I
Expected product life service	5 years
Caretaker must be able to physically and mentally operate the wheelchair safely, e.g. to drive, control and operate it	Yes
Temperature (storage and operation)	-10° C til 60° C
The surface temperature of the wheelchair can increase if exposed to high heat, e.g. strong sunlight over a longer period of time. The surface can also be cooled down in extreme cold	Yes
Is the wheelchair suitable as seat in vehicle?	No
If the basket under the seat is removed, will it destabilize the wheelchair?	No
Is the wheelchair (without user) suitable for land and air transport?	Yes
Produced in accordance with	EN 12182:2000
User guide is available in pdf form on website and can be scaled up to readable size for the visually impaired	Yes

For transporting bariatric users indoor at e.g. hospitals.

Model XXL-Rehab Transporter			
Seat width	Seat depth	Order nr	UDI DI
71 cm	56 cm	0182-071-000	5740000102803

