

Transportation of passenger in COBI MobiBed inside vehicle

This concerns the COBI / LASAL MobiBed , transportation of a person on a MobiBed and inside a moving vehicle.

Following applies:

1. COBI REHAB / LASAL CAN NOT RECOMMEND transportation of a person on a MobiBed inside a moving vehicle. Not even when the MobiBed is securely fastened during transport.
2. The MobiBed has not been subjected to crash testing. It is NOT possible to Crash Test a product such as the Mobibed according to the relevant standards.

According to standard: **E/ECE/324/Rev.1/Add.15/Rev.9** and **E/ECE/TRANS/505/Rev.1/Add.15/Rev.9** , the Crash Test can not be performed according to the safety purpose in the standard, due to seating angles and safety belt fixation positions.

According to standard: **ISO7176.19** , the Crash Test can not be performed according to the safety purpose in the standard, due to seating angles and safety belt fixation positions.

See also the Manual of the MobiBed product for correct handling.

For special cases ONLY: COBI REHAB recommends that a local professional evaluation of safety and risk related to the benefits, are made before people transportation purpose. The evaluator must assure that the amount of benefit clearly outweighs the amount of risk. COBI REHAB do NOT take any responsibility in such case.

COBI REHAB

September, 2020



Claus Nymark