



## WHEN NOTHING IS STRAIGHTFORWARD

### KAMILLE BACK SUPPORT FOR WHEELCHAIRS

The Kamille Back Support is unique because it can be adjusted asymmetrically. This makes it easier to address anatomical challenges such as scoliosis, increased lateral flexion in the trunk, or asymmetrical tissue, for example, due to overweight, hernias, or edema.

The Kamille Back Support ensures comfort, stability, and functionality for wheelchair users in need of a strong back support.

User weight up to 140 kg.



- Asymmetric adjustment of the backrest
- Adjustable side supports
- Positioning of the back can be adjusted in height and depth
- High recyclability, as it can be readjusted for a new user
- Safe and stable positioning solution
- Can be mounted with or without body supports
- Compatible with most standard wheelchairs
- Can be equipped with various types of neck supports

## WHEN NOTHING IS STRAIGHTFORWARD

### SPECIFICATIONS



Kamille Back Support is wheelchair back support with unique adjustment options. Three aluminium plates in the backrest and six positioning cushions, all of which can be adjusted individually.

The recyclability of the Kamille Back Support is very high. This is due to a combination of quality materials, the possibility of readjustment for a new user, and additional purchases.

The Kamille Back Support can be mounted on most electric and manual wheelchairs.

The height of the back is 500 mm + 200.5 mm depending on the placement of the back, while the width of the back is between 340 - 520 mm.

The cushion filling consists of viscoelastic foam and EPS beads. Cobi Rehab recommends purchasing a breathable and washable cover.

The grey, washable, and stretchable cover is coated with incontinence material on the back and moisture-absorbing PU-impregnated viscose on the outside.

The backrest is compatible with other neck supports.

The Kamille Back Support can be mounted by a consultant or a therapist on-site or by a technician at an assistive aid's depot.

**Contact Cobi Rehab at tel. +45 7025 2522 or at [cobi@cobi.dk](mailto:cobi@cobi.dk) for more information on item numbers, sizes, and installation.**



Scan the QR code and read more about the Kamille Back Support on our website.